

Sample Shipping Procedure and Pricing



Shipping Procedures

Monday through Thursday

- Place samples and ice pack in insulated envelope or bag prior to packing in box.
- Samples should be wrapped tightly and sealed to avoid samples being loose in the box and to ensure samples remain cool during shipping.

Friday

- Pack with extra ice packs to maintain samples over the weekend.

REGULAR SHIPPING METHOD

UPS, FedEx, DHL. US Postal

Use overnight or 2 day services only.
If shipping two day only ship Mon – Wed

Send to:

Animal Profiling International, Inc
6040 North Cutter Circle Suite 317
Portland, OR 97217
503-247-8066

1. Include client form with shipment
2. Sender will arrange payment

FREQUENT SHIPPING METHOD

UPS (discounted overnight or 2 day rates)

**API will provide preprinted shipping labels
Call for labels**

1. Include client form with shipment
2. Shipping will be billed on testing invoice

For questions contact: API Sales at (877) 278-1344 or API Lab at (503) 247-8066
Headquarters: 6040 N. Cutter Circle, Suite. 317, Portland, OR 97217 www.animalprofiling.com