

## SIMPLE STEPS TO TEST YOUR CATTLE

1. Visit our website or call API for sample collection and shipping protocols.
2. Collect a milk or blood sample according to API instructions.
3. Ship the samples to the API laboratory according to instructions.
4. Results will be reported directly to you and/or your veterinarian.

## ABOUT US

Animal Profiling International, Inc. is an American company dedicated to the improvement of animal health management through risk-assessment technologies.

API's industry experts provide the products and services that operators in the beef, dairy and food supply industries need to make their jobs more efficient and profitable.

The company has tested over 2.5 million head of cattle and dairy animals for BVD and currently services clients in over 40 U.S. states. In addition to BVD and pregnancy testing, the company also provides diagnostic services for the detection of BLV, Johne's Disease, Mastitis and Neospora.

[www.animalprofiling.com](http://www.animalprofiling.com)

# EarlyPreg28 *(blood pregnancy test)* EasyPreg *(milk pregnancy test)*

The **EarlyPreg28** and **EasyPreg** tests offer the unique opportunity to detect pregnancy earlier and with greater accuracy over traditional pregnancy testing methods while also reducing the risks of pregnancy loss due to palpation or ultrasound trauma. Additionally, the tests provide producers and operators further time flexibility, allowing them to test for pregnancy on their own schedules.

## EARLYPREG28

EarlyPreg28 is a blood-based pregnancy test designed to detect pregnancy in cattle and goats as early as 28 days post-insemination without invasive palpation.

The benefits to using the EarlyPreg28 test include:

- Accurate results 28 days post-breeding
- Convenience and versatility, allowing you to schedule pregnancy testing when it fits your schedule
- Non-invasive and embryo safe, reducing the risk of damaging a pregnancy and causing embryo loss

## EASYPREG

EasyPreg is an accurate milk pregnancy test designed to confirm or detect pregnancy in cattle and goats as early as 28 days post-insemination without invasive palpation.

The benefits to using the EasyPreg test include:

- The use of milk samples...no needle or special training required
- Schedule testing when it is right for you and the animals
- Confirm and dry-off checks on the dairy can easily be taken before milking
- Milk sample can be fresh, frozen or preserved

\* EarlyPreg28 is performed with the IDEXX Bovine Pregnancy Test.

\* EasyPreg is performed with the IDEXX Milk Pregnancy Test.

**FOR MORE INFORMATION, CONTACT API CUSTOMER SERVICE:  
877.278.1344 or [sales@animalprofiling.com](mailto:sales@animalprofiling.com)**